



Product Matrix



Normal Corneal Conditions	1st Choice	2nd Choice	3rd Choice
Spherical Error: Myopia and Hyperopia	D SimplifEyes 1Day Daily Disposable	H Duette	H SynergEyes A
Astigmatism up to 2.25D	H Duette	H SynergEyes A	GP SynergEyes GP II Bitoric
Astigmatism 2.5D or more	H Duette Enhanced Profile	H SynergEyes A Enhanced Profile	GP SynergEyes GP Bitoric
Residual Astigmatism > 1.0D	GP SynergEyes GP Front Toric	GP SynergEyes GP II Bitoric	
Presbyopia with Astigmatism up to 2.00D	H Duette Progressive Near or Distance Zone	GP SynergEyes GP II Progressive	GP SynergEyes GP Progressive
Presbyopia with Astigmatism 2.25D or more	H Duette Progressive Enhanced Profile Near or Distance Zone	GP SynergEyes GP Bitoric Progressive	

Irregular Corneal Conditions	1st Choice	2nd Choice	3rd Choice
Emerging Keratoconus	H UltraHealth	H UltraHealth FC	S SynergEyes VS
Moderate Keratoconus (inferior cone)	H UltraHealth FC	H UltraHealth	S SynergEyes VS
Moderate Keratoconus (central cone)	H UltraHealth	S SynergEyes VS	H SynergEyes KC
Advanced Keratoconus	S SynergEyes VS	H UltraHealth	H SynergEyes KC
Early Pellucid Marginal Degeneration	H UltraHealth FC	H UltraHealth	S SynergEyes VS
Moderate Pellucid Marginal Degeneration	S SynergEyes VS	H UltraHealth FC	H UltraHealth
Advanced Pellucid Marginal Degeneration	S SynergEyes VS	H UltraHealth	
Keratoglobus	S SynergEyes VS	S SynergEyes VS	
Post Keratoplasty - Oblate Graft	S SynergEyes VS	H UltraHealth FC	
Post Keratoplasty - Prolate Graft	S SynergEyes VS	H UltraHealth	
Post Surgical - Rk	H UltraHealth FC	H UltraHealth	S SynergEyes VS
Post Surgical Ectasia (mild to moderate)	H UltraHealth FC	H UltraHealth	S SynergEyes VS
Post Surgical Ectasia (severe)	S SynergEyes VS	H UltraHealth	

Legend:

D Daily Disposable Soft*	H Hybrid Lenses Normal Cornea	S Scleral Lenses*	GP Corneal Gas Permeable Lenses*	H Hybrid Lenses Irregular Cornea
---------------------------------	--------------------------------------	--------------------------	---	---

*Daily Disposable Lenses and Gas Permeable Lenses available in US only.